

LÖDIGE: THE NEXT LEVEL IN COATING TECHNOLOGY

- 25% faster processes with outstanding coating uniformity
- 10% to 100% fill levels without drum or baffle changes
- Optimised air distribution = less spray drying, faster cleaning, reduced downtime
- Contactless sampling
- Turn-Key projects
- Fully customized

ALWAYS THE RIGHT MIX

LÖDIGE COATERS SUCCESSFUL SINCE 1980



Whether you're focused on scale-up, mini-tablet capability, or reducing production costs, the LC Coater delivers measurable improvements from lab to full-scale manufacturing.



Experience of more than 70 years in mixing were transferred in new baffle design for perfect mixing results.



Distance to tablet bed adjustable from outside (manually or fully automatic) without disturbing the coating process.



Automatic Sampling using a Pneumatic Vacuum System. No direct contact with the Product. Increased operator safety!

The following products can be processed with a suitable coating drum design:

- Tablets of all shapes
- Mini-tablets
- Sugar coated tablets
- Capsules
- Pellets



**BOOK A
TEAMS MEETING
OR ARRANGE A
TRIAL AT THE
LÖDIGE TEST
CENTRE**



Torsten Meinel-Dirumdam
Meinel-Dirumdam@loedige.de
Phone +49 5251 309278



Matthias Niggemann
Niggemann@loedige.de
Phone +49 5251 309312

www.loedige.de/en/process/coating